
TECHNICAL DATA SHEET

Device picture			
Product:	AIDAL	Reference Number:	437000
GMDN code:	<p>36777 - An assembly of devices used to evacuate gas, fluid, tissue, or other foreign materials from a body cavity or lumen by means of suction. It generally consists of a motorized suction pump, tubing, plastic/glass collection container(s), a vacuum gauge, a vacuum control knob, an overflow trap, a moisture filter, and possibly a microbial filter. The pump creates a vacuum in the suction tubing, which is inserted into the body for the removal of materials into the collection container. This system can be used in a wide variety of settings within healthcare facilities.</p>		
Medical device classification:	Class IIb		
CE marking:	Certificate n°679/MDD in accordance with the requirements of the European Directive 93/42/EEC and following modifications – Annex II		
Manufacturer:	Air Liquide Medical Systems S.r.l. Via dei Prati, 62 25073 Bovezzo (BS) – ITALIA		
Manufacturer certifications:	UNI CEI ISO 13485:2012 UNI EN ISO 14001:2015 BS OHSAS 18001:2007		
Product normative references:	UNI EN ISO 10079-1:2009		
Description:	<p>AIDAL is a surgical suction pump with light switch, vacuum gauge and flux regulator.</p> <p>It can be used with a reusable or single use jar of 1000ml.</p> <p>The aspirator can be operated on a rechargeable 12V DC battery or a standard mains power supply.</p> <p>The package contains:</p>		

	<ul style="list-style-type: none"> • The aspirator and rechargeable battery • A power supplier adapter and cord • A 1000 ml graduated collection jar • Two antibacterial filters • Tube kit (2 tubes of varying length) • A reducing connector to connect a sterile disposable cannula and the connection tube • Operating instructions and warranty • Aspirator carrying bag with strap • Connector with valve for suction control
Intended use:	The AIDAL surgical aspirator is a portable device intended for the suction of fluids and other substances from the human body. It can be used in the pharyngeal area and with tracheotomised patients. Vacuum pressure and flow can be adjusted. The AIDAL aspirator is not intended for thoracic drainage.

Materials:	COMPONENTS	MATERIAL
	External body	ABS V0 + Polypropylene
Jar	Polycarbonate	
Battery	Pb	
Tubings	PVC	
External box	Cardbord	
User manual	Paper	

Technical specifications:	<p>Voltage: 12V DC Maximum absorption: 4 A Battery: 12V DC 2.3Ah Pb Battery powered operating time: 40 minutes (depending on use) Switching adapter with cord: 100-240V 50/60Hz - 15V DC 4A Functioning: continuous Weight with jar: 3.8 kg Capacity of the standard jar: 1000 ml Dimensions: 350 x 140 x 280 (H) mm Maximum flow rate at the pump: 30±4 lpm Maximum vacuum: -80 kPa ± 8% Protection level IP: IP34 Electrical class II Part applied BF</p> <p>Flow-Vacuum diagram:</p> 
----------------------------------	---

	Measured with Fluke VT plus (error \pm 4 lpm)
Cleaning, disinfection and sterilization:	<p>CLEANING AND DISINFECTION: The aspirator can be cleaned and disinfected using a sponge or soft cloth moistened with a cold disinfectant solution (e.g. a sodium hypochlorite solution with a concentration no greater than 2%), as directed by the manufacturer of the product.</p> <p>STERILIZATION: Place the parts in an autoclave for a sterilization cycle with a steam temperature of 121°C for at least 15 minutes (preferably not exceeding 30 minutes) taking care to position the graduated jar upside down (with its bottom upward). Steam sterilization may be performed for a maximum of 15 cycles without adversely affecting the characteristics of the component materials used and therefore functionality. However, a visual inspection of all the components is recommended before each assembly following sterilization. Sterilization can also be carried out cold by ethylene oxide or gamma rays following the instructions supplied with respective equipment.</p>
User guide:	The user manual in each package English, Italian, French, German, Spanish, Portuguese, Polish, Dutch, Swedish, Finnish, Greek, Slovak, Japanese and Korean, Czech.
Storage conditions:	<ul style="list-style-type: none"> • Temperature -15 to +40 °C • Relative humidity 0 to 93% • Pressure 700 to 1060 hPa
Traceability:	Serial number
Packaging:	Single packaging
Date of revision:	Rev.D 10/2018


Air Liquide Medical Systems S.r.l.
Technical Department

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Air Liquide Medical Systems S.r.l.
Bureau technique